



PRO BASE®

CATALOGUE 2025



Prodenta3D –

a technological center that has been successfully operating for more than a decade. We started our activities by manufacturing custom metal structures for specific cases, and today we proudly present our long-developed and refined product – **PROBASE® titanium abutments**. This is an advanced prosthetic solution characterized by highly flexible application possibilities and adapted for digital workflows.

Our accumulated experience, continuous innovations, modern equipment, and a team of highly qualified specialists allow us to ensure that **PROBASE®** titanium abutments are not only innovative but also meet the strictest European safety standards. Striving for the highest quality, we offer not only advanced **PROBASE®** abutments but also comprehensive collaboration to help implement all prosthetic solutions.



Certificates – a guarantee of quality and safety

The **Prodenta3D** production unit meets the strictest quality and safety standards, ensuring top-quality manufacturing and reliable products. Our management system is certified according to the **ISO 13485 standard**, which regulates the manufacturing processes of medical devices. This certification ensures that every stage of production and management is strictly controlled, creating a safe environment and maintaining the highest quality standards.

Our **PROBASE®** titanium abutments comply with the latest Medical Device Regulation (MDR) requirements. This certification allows our products to be marked with the **CE mark**, confirming that they meet all European safety and quality requirements for medical devices.

The obtained certificates not only prove our commitment to quality but also confirm that **Prodenta3D** is a recognized European manufacturer of medical devices, meeting the highest international manufacturing standards.

Do we have guarantees? Yes! Not only for abutments but also for implants!

Prodenta3D guarantees the highest quality and durability – that's why we offer a lifetime warranty on all **PROBASE®** abutments. This means that by choosing our products, you can be completely confident in their reliability and longevity.

However, our commitments do not end there. We also provide a warranty for the dental implants on which prosthetics were

performed using **PROBASE®** atramas. Ši garantija atitinka gamintojo nustatytas sąlygas.

Pasirinkdami **Prodenta3D** abutments. This warranty complies with the conditions set by the manufacturer.

By choosing Prodenta3D, you can rest assured – we take full responsibility for our products and ensure that they will last a lifetime.



**Lifetime warranty
for Probase®
abutments**



**Third-party
warranty for
dental implants**

Advantages of Probase® Abutments



Probase® Multi-Unit products

Multi-unit abutment

5,5 mm

4,5 mm

Two working heights (compatible with CAD/CAM libraries)

Libraries created for digital workflows

Designed for direct zirconia screw-retention on Multi-unit abutments

The large cone provides a strong grip

We provide complete libraries and manufacturing strategies

R-Type screw

Universal scanning body

Large scanning area

Libraries created for various solutions. With direct connections or through abutments.

Titanium connection ensures stability

Compatible with most Multi-Unit systems

CONTENTS	
Implant level	7
Connical ANK Implant <i>Compatible with Ankylos® Implant</i>	7
Connical EV M <i>Compatible with Astra Tech® EV M</i>	8
Connical EV S <i>Compatible with Astra Tech® EV S</i>	9
Connical AR <i>Compatible with Megagen® Anyridge</i>	10
Connical MGG AO <i>Compatible with Megagen® AO</i>	11
Connical MGG Octa <i>Compatible with Megagen® Octa</i>	11
Connical GM Ø3.8 <i>Compatible with Neodent® GM Ø3.8</i>	12
Connical GM Ø4.5 <i>Compatible with Neodent® GM Ø4.5</i>	13
Connical NP <i>Compatible with Nobel® NP</i>	15
Connical RP <i>Compatible with Nobel® RP</i>	16
Connical Mini <i>Compatible with Osstem® TS mini</i>	18
Connical Regular TS <i>Compatible with Osstem® TS regular</i>	19
Connical NC <i>Compatible with Straumann® NC Bone level</i>	21
Connical RB <i>Compatible with Straumann® RB Bone level</i>	22
Connical RC <i>Compatible with Straumann® RC Bone level</i>	23
Octagon RN <i>Compatible with Straumann® RN Tissue level</i>	25
Multi-unit level	26
Multi ANK <i>Ankylos® MU 4.2 scanbody</i>	26
Multi MED <i>Compatible with Medentika® Multi-unit</i>	27
Multi MGG <i>Compatible with Megagen® Multi-unit</i>	28
Multi NEO <i>Compatible with Neodent® Multi-unit</i>	29
Multi NOB <i>Compatible with Nobel® Multi-unit ø4.8</i>	30
Multi OSS <i>Compatible with Osstem® Multi-unit</i>	31
Multi STR <i>Compatible with Straumann® Multi-unit ø4.6</i>	32
Universal Multi-unit line	34

Implant level



Connical ANK Implant

Compatible with Ankylos® Implant



Probase® screw			
REF. Code	Screwdriver	Quantity in the package	Torque force
17.01.SC00X1	Hex 1.0	1 pcs.	35Ncm
17.01.SC00X5	Hex 1.0	5 pcs.	35Ncm

Connical EV M

Compatible with Astra Tech® EV M

Probase® abutment for crown

- Titanium
- Torque force 25Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm
Rotation	Engaging		
Screw	19.01.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	19.01.AB01X3	19.01.AB02X3	19.01.AB03X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
19.01.SC00X1	Hex 1.26	1 pcs.	25Ncm
19.01.SC00X5	Hex 1.26	5 pcs.	25Ncm

Connical EV S

Compatible with Astra Tech® EV S

Probase® abutment for crown

- Titanium
- Torque force 25Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm
Rotation	Engaging		
Screw	18.01.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	18.01.AB01X3	18.01.AB02X3	18.01.AB03X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
18.01.SC00X1	Hex 1.26	1 pcs.	25Ncm
18.01.SC00X5	Hex 1.26	5 pcs.	25Ncm

Connical AR

Compatible with Megagen® Anyridge

Probase® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm	4mm
Rotation	Engaging			
Screw	18.01.SC00			
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	04.01.AB01X3	04.01.AB02X3	04.01.AB03X3	04.01.AB04X3

Probase® abutment for bridge

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm	4mm
Rotation	Non Engaging			
Screw	04.01.SC00			
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	04.01.AB11X3	04.01.AB12X3	04.01.AB13X3	04.01.AB14X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
04.01.SC00X1	Hex 1.2	1 pcs.	35Ncm
04.01.SC00X5	Hex 1.2	5 pcs.	35Ncm

Connical MGG AO

Compatible with Megagen® AO

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
04.04.SC00X1	Hex 1.2	1 pcs.	35Ncm
04.04.SC00X5	Hex 1.2	5 pcs.	35Ncm

Connical MGG Octa

Compatible with Megagen® Octa

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
04.02.SC00X1	Hex 1.2	1 pcs.	35Ncm
04.02.SC00X5	Hex 1.2	5 pcs.	35Ncm

Connical GM Ø3.8

Compatible with Neodent® GM Ø3.8

Probased® abutment for crown

- Titanium
- Torque force 20Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	0.8mm	1.5mm	2.5mm
Rotation	Engaging		
Screw	03.01.SC01		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	03.01.AB01X3	03.01.AB02X3	03.01.AB03X3

Probased® abutment for bridge

- Titanium
- Torque force 20Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	0.8mm	1.5mm	2.5mm
Rotation	Non Engaging		
Screw	03.01.SC01		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	03.01.AB11X3	03.01.AB12X3	03.01.AB13X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
03.01.SC01X1	Unigrip	1 pcs.	20Ncm
03.01.SC01X5	Unigrip	5 pcs.	20Ncm

Connical GM Ø4.5

Compatible with Neodent® GM Ø4.5

Probased® abutment for crown

- Titanium
- Torque force 20Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	0.8mm	1.5mm	2.5mm
Rotation	Engaging		
Screw	03.01.SC01		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	03.01.AB06X3	03.01.AB07X3	03.01.AB08X3

Probased® abutment for bridge

- Titanium
- Torque force 20Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	0.8mm	1.5mm	2.5mm
Rotation	Non Engaging		
Screw	03.01.SC01		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	03.01.AB16X3	03.01.AB17X3	03.01.AB18X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
03.01.SC01X1	Unigrip	1 pcs.	20Ncm
03.01.SC01X5	Unigrip	5 pcs.	20Ncm

Connical NP

Compatible with Nobel® NP

Probase® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm
Rotation	Engaging		
Screw	02.02.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	02.02.AB01X3	02.02.AB02X3	02.02.AB03X3

Probase® abutment for bridge

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm
Rotation	Non Engaging		
Screw	02.02.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	02.02.AB04X3	02.02.AB05X3	02.02.AB06X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
02.02.SC00X1	Unigrip	1 pcs.	35Ncm
02.02.SC00X5	Unigrip	5 pcs.	35Ncm

Connical RP

Compatible with Nobel® RP

Probase® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm
Rotation	Engaging		
Screw	02.03.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	02.03.AB01X3	02.03.AB02X3	02.03.AB03X3

Probase® abutment for bridge

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1mm	2mm	3mm
Rotation	Non Engaging		
Screw	02.03.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	02.03.AB04X3	02.03.AB05X3	02.03.AB06X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
02.03.SC00X1	Unigrip	1 pcs.	35Ncm
02.03.SC00X5	Unigrip	5 pcs.	35Ncm

Connical Mini

Compatible with Osstem® TS mini

Probase® abutment for crown

- Titanium
- Torque force 20Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Engaging		
Screw	07.03.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	07.03.AB01X3	07.03.AB02X3	07.03.AB03X3

Probase® abutment for bridge

- Titanium
- Torque force 20Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Non Engaging		
Screw	07.03.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	07.03.AB04X3	07.03.AB05X3	07.03.AB06X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
07.03.SC00X1	Hex 1.2	1 pcs.	20Ncm
07.03.SC00X5	Hex 1.2	5 pcs.	20Ncm

Connical Regular TS

Compatible with Osstem® TS regular

Probased® abutment for crown

- Titanium
- Torque force 30Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Engaging		
Screw	07.02.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	07.02.AB01X3	07.02.AB02X3	07.02.AB03X3

Probased® abutment for bridge

- Titanium
- Torque force 30Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Non Engaging		
Screw	07.02.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	07.02.AB05X3	07.02.AB06X3	07.02.AB07X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
03.01.SC01X1	Unigrip	1 pcs.	20Ncm
03.01.SC01X5	Unigrip	5 pcs.	20Ncm

Connical NC

Compatible with Straumann® NC Bone level

Probase® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Engaging		
Screw	01.05.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	01.04.AB01X3	01.04.AB02X3	01.04.AB03X3

Probase® abutment for bridge

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Non Engaging		
Screw	01.05.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	01.04.AB11X3	01.04.AB12X3	01.04.AB13X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
01.05.SC00X1	Torx T6	1 pcs.	35Ncm
01.05.SC00X5	Torx T6	5 pcs.	35Ncm

Connical RB

Compatible with Straumann® RB Bone level

Probase® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.3mm	2.5mm	3.5mm
Rotation	Engaging		
Screw	01.01.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	01.01.AB04X3	01.01.AB02X3	01.01.AB03X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
01.01.SC00X1	Torx T6	1 pcs.	35Ncm
01.01.SC00X5	Torx T6	5 pcs.	35Ncm

Connical RC

Compatible with Straumann® RC Bone level

Probased® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Engaging		
Screw	01.05.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	01.05.AB01X3	01.05.AB02X3	01.05.AB03X3

Probased® abutment for bridge

- Titanium
- Torque force 30Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	1.0 mm	2.0 mm	3.0 mm
Rotation	Non Engaging		
Screw	01.05.SC00		
Abutment heights	5.5/4.3/3.2	5.5/4.3/3.2	5.5/4.3/3.2
REF. code	01.05.AB11X3	01.05.AB12X3	01.05.AB13X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
01.05.SC00X1	Torx T6	1 pcs.	35Ncm
01.05.SC00X5	Torx T6	5 pcs.	35Ncm

Octagon RN

Compatible with Straumann® RN Tissue level

Probase® abutment for crown

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	-
Rotation	Engaging
Screw	01.07.SC00
Abutment heights	6/4.85/3.7
REF. code	01.04.AB01X3

Probase® abutment for bridge

- Titanium
- Torque force 35Ncm
- **Screw included**
- 3 Abutment heights



Gingiva height	-
Rotation	Non Engaging
Screw	01.07.SC00
Abutment heights	6/4.85/3.7
REF. code	01.04.AB11X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
01.07.SC00X1	Torx T6	1 pcs.	35Ncm
01.07.SC00X5	Torx T6	5 pcs.	35Ncm

Multi-unit level



Multi ANK

Ankylos® MU 4.2 scanbody

Probase® scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw	17.02.SC00
REF. code	17.02.SB03

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
17.02.SC00X1	Hex 1.0	1 pcs.	15Ncm
17.02.SC00X5	Hex 1.0	5 pcs.	15Ncm

Multi MED

Compatible with Medentika®Multi-unit

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
17.02.SC00X1	Hex 1.0	1 pcs.	15Ncm
17.02.SC00X5	Hex 1.0	5 pcs.	15Ncm

Probased® Universal MU Ø4.6-Ø4.8 scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw

01.03.SC00

REF. code

23.01.SB03

**Probased® R-Type Screw**

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Multi MGG

Compatible with Megagen® Multi-unit

Probased® abutment for bridge

- Titanium
- Torque force 15Ncm
- **Screw included**
- 2 Abutment heights



Screw

04.03.SC00

Abutment heights

5.65/4.5

REF. code

04.03.AB01X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
04.03.SC00X1	Hex 1.2	1 pcs.	15Ncm
04.03.SC00X5	Hex 1.2	5 pcs.	15Ncm

Probased® Universal MU Ø4.6-Ø4.8 scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw

01.03.SC00

REF. code

23.01.SB03

**Probased® R-Type Screw**

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Multi NEO

Compatible with Neodent® Multi-unit

Probase® abutment for bridge

- Titanium
- Torque force 15Ncm
- **Screw included**
- 2 Abutment heights



Screw	03.02.SC00
Abutment heights	5.65/4.5
REF. code	03.02.AB01X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
03.02.SC00X1	Unigrip	1 pcs.	15Ncm
03.02.SC00X5	Unigrip	5 pcs.	15Ncm

Probase® Universal MU Ø4.6-Ø4.8 scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw	01.03.SC00
REF. code	23.01.SB03



Probase® R-Type Screw

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Multi NOB

Compatible with Nobel® Multi-unit Ø4.8

Probase® abutment for bridge

- Titanium
- Torque force 15Ncm
- **Screw included**
- 2 Abutment heights



Screw	02.01.SC01/ 02.01.SC00
Abutment heights	5.65/4.5
REF. code	02.01.AB01X3



Probase® Screw

REF. Code	Screwdriver	Quantity in the package	Torque force
02.01.SC01X1	Unigrip	1 pcs.	15Ncm
02.01.SC01X5	Unigrip	5 pcs.	15Ncm



REF. Code	Screwdriver	Quantity in the package	Torque force
02.01.SC00X1	Omnigrip	1 pcs.	15Ncm
02.01.SC00X5	Omnigrip	5 pcs.	15Ncm

Probase® Universal MU Ø4.6-Ø4.8 scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw	01.03.SC00
REF. code	23.01.SB03

**Probased® R-Type Screw**

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Multi OSS

Compatible with Osstem® Multi-unit

Probased® abutment for bridge

- Titanium
- Torque force 20Ncm
- **Screw included**
- 2 Abutment heights



Screw	07.01.SC00
Abutment heights	5.65/4.5
REF. code	07.01.AB01X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
07.01.SC00X1	Hex 1.2	1 pcs.	20Ncm
07.01.SC00X5	Hex 1.2	5 pcs.	20Ncm

Probased® Universal MU Ø4.6-Ø4.8 scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw	01.03.SC00
REF. code	23.01.SB03

**Probased® R-Type Screw**

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Multi STR

Compatible with Straumann® Multi-unit Ø4.6

Probased® abutment for bridge

- Titanium
- Torque force 15Ncm
- **Screw included**
- 2 Abutment heights



Screw	01.03.SC00
Abutment heights	5.65/4.5
REF. code	01.03.AB01X3

**Probased® Screw**

REF. Code	Screwdriver	Quantity in the package	Torque force
01.03.SC00X1	Torx T6	1 pcs.	15Ncm
01.03.SC00X5	Torx T6	5 pcs.	15Ncm

Probased® Universal MU Ø4.6-Ø4.8 scanbody

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw	01.03.SC00
REF. code	23.01.SB03




Probaser® R-Type Screw

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Universal
Multi-
unit line



Probase® Universal MU Ø4.6-Ø4.8 scanbody

Compatible with Medentika® Multi-unit, Nobel® Multi-unit Ø4.8, Neobiotech® Multi-unit, Osstem® Multi-unit, Neodent® Multi-unit, Straumann® Multi-unit Ø4.6, Megagen® Multi-unit systems

- Titanium + PEEK
- Torque force 15Ncm
- **Screw included**



Screw	01.03.SC00
REF. code	23.01.SB03



Probase® R-Type Screw

Designed for direct zirconia/PMMA screwing

REF. Code	Screwdriver	Quantity in the package	Torque force
23.01.SC00X1	Torx T6	1 pcs.	15Ncm
23.01.SC00X5	Torx T6	5 pcs.	15Ncm

Ankylos® is a registered trademark owned by Dentsply Sirona Implants.
Astra Tech® is a registered trademark owned by Dentsply IH AB.
Megagen® is a registered trademark owned by MEGAGEN IMPLANT CO. LTD.
Neodent® is a registered trademark owned by JJGC Indústria e Comércio de Materiais Dentários S.A.
Nobel® is a registered trademark owned by Nobel Biocare AB.
Osstem® is a registered trademark owned by OSSTEM IMPLANT CO. LTD.
Straumann® is a registered trademark owned by Straumann Holding AG, Switzerland.

CONTACTS

orders@prodenta.it

+370 6 568 0590

www.prodenta3d.com



prodenta 3D